Physical Rehabilitation Lower limbs Functional Biofeedback Analytical Eccentric Isokinetics Evaluation Follow-up of the patient practice Exercise control Upper limbs Therapeutic strategies Concentric Passive mobilisation

# A new generation of aid: "all-in-one"!

## A product range designed by and for physiotherapists to meet their needs and those of their patients!

The second generation **Kinevolution**, well-designed and compact, has been developed with the aid of physiotherapists-massage therapists, the majority of whom are in private practice. Thanks to this range of assistants, the physiotherapist sets up **multiple therapeutic strategies** for the purpose of prevention, maintenance, rehabilitation, recovery and evaluation of **all of his patient practice**!

The therapist simply starts, on a large dimension **touch screen**, a session of rehabilitation with exercises for **control of exertion** (re-creation of mechanical stresses of manual resistance), with analytical or **functional isokinetic exercise**, with **passive mobilisation** or **isometric exercise**.

The practitioner and his patient visualise, **in real time** on this new interface, **biofeedback** of exertion performed compared to the specification that the apparatus quantifies with **precision**. Patient follow-up is facilitated by using the medical dossier, which can be consulted at all times on the equipment, and printable results can shared with the medical team.

**Kin**evolution meets the needs of physiotherapists in Hospital Centres, in Functional Rehabilitation Centres, in Health Care and Rehabilitation institutions, so that they can adapt to the current progression of the profession without limitations.

## An assistant available throughout the healthcare course!

**Kinevolution** is a multi-therapeutic equipment which offers physiotherapists the **widest range of therapeutic uses** for the patient's health care course in orthopaedics, rheumatology, traumatology or in neurology.

The therapist also benefits from measured performance, which justifies the results of each session. Thus, he/she performs objective follow-up of his/her patient practice in the centre.

### A gain in time and energy at each session!

The therapist is **freed** from the strain of performing manual resistance and can **focus on the patient** whose **cognitive functions** are stimulated by real time biofeedback for more effective **coaching!** 

In less than **40 seconds** it is possible to schedule a session comprised of several exercises that the physiotherapist readily initiates with **fingertip control** on the touch screen!



Visit our website:

#### www.kinevolution.com

Or contact us: contact@kinevolution.com +334 66 549 568 +336 68 203 544

## Certification: a guarantee of quality and safety!

То satisfy the needs of physiotherapists and massage therapists in the best possible manner in treatment centres. Kinevolution has CE certification as a class II medical device, thus justifying quality, safety and regulatory compliance of our equipment.

Similarly, our design and production operations, performed exclusively in France, meet the requirements of the **ISO 13485** standard.

## Innovation dedicated to physical rehabilitation